H&O EQUIPMENTS  
Reinventing Cryosurgery

User Manual


Cryo**Probe**™ B+

## Table of contents

User manual	3-17
Notice explicative	18-32
Instrucciones de uso	33-47

## User manual

**SN** The serial number is engraved in the cavity.  
Please mention your instrument's serial number  
in all communication with H&O Equipments.

This user manual is applicable  
for the following products and accessories.

## Instruments

- *Expected lifetime : minimum 10 years in normal use.*

With orange/red, white/red, green/red applicators

**Cartridges** - *Expired date on the cartridges packaging (shelf life three years).*

23,5g N<sub>2</sub>O cartridges: 1 carton of 12 boxes  
each box containing 6 cartridges

## Applicators

Orange/red dot applicator for 1-4 mm applications  
White/red dot applicator for 2-6 mm applications  
Green/red dot applicator for 4-8 mm applications

## Accessories

Tweezers for filter/opener placement  
Locking cap for applicators

**REF**

## CryoProbe B+

S-HO-CBPLUSPR-05

S-HO-NO23-12-S06-PR

S-HO-CCX0-MA-027-PR  
S-HO-CCX0-MA-005-PR  
S-HO-CCX0-MA-006-PR

S-HO-AXXX-FO-001  
S-HO-APPL-CAP

The above list of products is correct at the time of printing. Some articles may no longer be available when reading this document. Some articles may not be available in certain countries.

EC or Canadian Declaration of Conformity is available and can be obtained upon simple request.

CryoProbe™ is the registered trademark of H&O Equipments nv/sa and is distributed under that name in the US and Canada.

CryoProbe™ is manufactured by H&O Equipments nv/sa.

FDA certificate is available and can be obtained upon simple request.

CryoProbe instruments are intended for the controlled destruction of unwanted tissue by application of extreme cold using liquefied N<sub>2</sub>O. Pathological indications are benign skin lesions such as: Acrochordon, actinic keratosis and cheilitis, capillary haemangioma, cutaneous leishmaniasis, discoid lupus erythematosus, filiform, flat, vulgar, plantar and anogenital warts, rosacea, lichen sclerosus, granuloma annulare, molluscum contagiosum, seborrheic keratosis, viral lesions by dispensing an extremely precise flow of nitrous oxide applied directly on the lesion.

All mentioned products are manufactured in the European Community by:

 H&O Equipments nv/sa, Rue des Journaliers 1,  
7822 Ghislenghien, Belgium

This manual's version: 200226CBPR

Telephone: +32 68 26 86 00

Fax: +32 68 55 28 48


E-mail: [info@ho-equipments.com](mailto:info@ho-equipments.com)

Website: [www.ho-equipments.com](http://www.ho-equipments.com)

The most recent version of this user manual and other language versions can be downloaded on our website.

Attention! Reference is made to the liability waiver at the back of this manual.

## Used symbols

-  Manufacturer
-  Caution, consult accompanying documents
-  Consult the user's guide
-  Do not expose the cartridges to temperatures above 50 °C/122 °F.

Keep away from sunlight

Catalogue number

Serial number

Batch number

Use-by date

Notified body nr 1639  
SGS BELGIUM

## IMPORTANT NOTICE

- When used correctly, the device does not come into contact with the lesion nor the skin.
- When used normally, there is no need to perform any maintenance, cleaning or disinfection.
- Should the applicator touch the patient in case of misuse, the applicator can be cleaned (autoclaved: 134 °C/273 °F - 4min or 121 °C/250 °F – 20 min) and disinfected (use a non-corrosive disinfecting solution, follow the disinfecting solution manufacturer's recommendations with regards to the concentration and time of application).
- Read the security instructions and remove the information tag from the instrument body.
- Always gently manipulate all components of CryoProbe instruments. Never use more than gentle hand force. Never use pliers or other tools to manipulate the instrument.
- Please ensure to only use the filter/opener supplied with the N<sub>2</sub>O cartridges. Any other filter/opener would damage the instrument.
- Incorrect placement of the filter/opener device into the loading cavity will cause irreversible damage to the instrument.
- Residual foreign matters in the loading cavity while placing a new filter/opener can cause instrument failure.
- When loaded, the instrument is under high pressure. Any change to the construction of the instrument may cause a potential risk.
- Never put your hand on the opening while removing the cartridge!
- All external parts can be wiped with a cloth soaked in any noncorrosive disinfecting solution. Never immerse the entire instrument in a disinfecting solution. Only the applicators can be autoclaved.
- The applicators and the cartridges must be kept away from dust and high temperatures. Leave the cap on when not in use.
- The CryoProbe should be handled with care. When storing the instrument with a loaded cartridge in its case, position the CryoProbe in such a way that the trigger switch cannot be depressed when closing the lid of the case. This would cause undesired gas flow.
- CryoProbe instruments are intended for professional use only.
- The N<sub>2</sub>O cartridges contain both a liquid and a gaseous phase in balance. The proportion of each phase depends on the room temperature which affects the output flow. The best performances are achieved when the cartridge temperature is between 19 and 22°C (66 and 72°F).
- Empty cartridges should be discarded as metal scrap.

## **Contraindications :**

- Unstable diabetes
- Skin conditions, e.g. skin tumors, esanthema, open wounds, solar hyperkeratosis
- Unexplained, suspicious liver spots and moles
- Cancer tissues and malignant tumors

## **The contraindications related to temporary conditions :**

- Infections accompanied by fever
- Acute chemotherapy or radiotherapy from four weeks before beginning the therapy to four weeks after finishing the therapy
- Pregnancy or breastfeeding
- Cold intolerance
- Vascular insufficiency

## **Reduce possible side-effects**

Although cryotherapy is a relatively low-risk procedure, some side effects may occur as a result of the treatment. They include:

- **Permanent changes in pigmentation:** Both hypopigmentation (lightening of the skin) and hyperpigmentation (darkening of the skin) may occur after cryotherapy. Both generally last a few months, but can be longer lasting. Avoid freezing the basal cell layer where melanocytes (pigment producing cells) are located.
- **Sensory impairment:** Though rare, damage to nerves is possible, particularly in areas where they lie closer to the surface of the skin, such as the fingers, the wrist, and the area behind the ear. Reports suggest this will disappear within several months.
- **Spattering of the cryogen during spraying, when the end nozzle freezes:** The innovation of CryoProbe is the direct application of nitrous oxide under high pressure (50-55 bar). This high pressure jet may cause minor shards of frozen humidity in the air blown away in a circle of approximately 30cm of diameter. They will thaw the moment they would eventually touch healthy skin.
- **Hair loss:** Hair follicles are easily damaged by cryosurgery and permanent hair loss is not uncommon.
- **Pigmented spots. It is common use to not treat pigmented spot in sunny seasons:** Solar lentigines (SL) are flat brown sharply demarcated pigmented lesions that commonly appear as a sign of aging and photodamage on sun-exposed skin.

- **Post-inflammatory hyperpigmentation (PIH)** is a frequent concern when treating SL in darker skin types.

Cryotherapy using CryoProbe is a simple and practical therapy for SL considering the following instructions:

- 1- Select your candidate: A typical candidate for treatment has multiple, brown, well defined, small and homogenous color lesions with no history of significant post-inflammatory hyperpigmentation.
- 2- Perform a spot test by treating one lesion, and evaluate the result in three weeks.
- 3- Apply sunscreen and use whitening cream for darker skin types for 3 weeks before and after treatment.

### **Further side-effects**

- Edemas
- Bleeding
- Pain and syncope
- Insufflation of subcutaneous tissue
- Retraction caused by freezing next to orifices
- Cartilage defect during treatment of the nose or the ear
- Infection
- Conjunctival ectopy due to treatments of the eyelid
- Milia
- Depressed or atrophic scars

## Step 1 | Empty the loading cavity

Before loading a new cartridge, release any unused gas from the used cartridge by unscrewing the applicator tip followed by pressing the button. Discard the used cartridge. Please ensure that the used filter does not remain inside the cavity.

### **Important note:**

In the case that you would not follow the gas release procedure as described before, the discharge of the gas will pass through the small opening under the on/off button. Extreme attention should be paid to avoid directing the jet of escaping gas to the user's hands, face or any other unprotected body parts. Never cover this opening! Never put your hand on this opening while removing the cartridge! This action must be done in a well ventilated room and by preference at an open window.

Empty loading cavity


Blistered cartridge and assembled filter/opener device

## Step 2 | Load a new cartridge

The cartridges are packed in peel-pack blisters with an assembled filter and cartridge opener included. The filters are used for protection of the applicator. The opener punctures the cartridge.

CryoProbe B+ can be used with 23,5g cartridges only. Ensure that the loading cavity is empty of all residue.

Open the peel-pack. Take out the filter/opener with the tweezers. Hold the CryoProbe's tip end pointed downward. Place the filter in the instrument and make sure that it is seated flat at the bottom of the loading cavity.

## Step 3 | Puncture the unused cartridge

Gently slide the cartridge into the CryoProbe's loading cavity and start screwing the cartridge turning it clockwise in a quick and uninterrupted movement until hand-tight.

Installation of the assembled filter/opener device

#### Step 4 | Use your CryoProbe B+

Remove the applicator cap. Gently squeeze the trigger switch to start the working. Releasing the trigger will close off the flow of cryogenic liquid. You may observe that during 1 or 2 seconds after releasing the trigger, the remaining gas in the applicator is discarded.

Applicators with various flow rates are available. Ask your distributor or check on [www.ho-equipments.com](http://www.ho-equipments.com).

	23,5g	Lesion size	Power
● Orange/Red	375-400 s	1-4 mm	▬
○ White/Red	200-250 s	2-6 mm	▬▬
● Green/Red	150-170 s	4-8 mm	▬▬▬

## Position the CryoProbe for treatment

The instrument must be held vertically with the tip directed downward. An angle of less than 45° must be avoided.

## Effective application

Before starting treatment take into account age, location, size and number of lesions, skin type, degree of tanning and season.

The applicator tip (from where the liquefied N<sub>2</sub>O flows) should be held 0,5 to 4 mm from the lesion. Further than 5 mm from the applicator tip the liquefied N<sub>2</sub>O will have expanded into gas phase which is not effective for achieving therapeutic results.


By approaching the lesion with the gas flowing you will observe the focal point of the jet on the tissue. The closer you approach the lesion, the more freezing power you apply to the tissue. It is recommended to keep a greater distance for the smaller and more delicate procedures.

The amount of liquid N<sub>2</sub>O that you wish to apply can be controlled by moving the applicator quickly towards and away from the tissue.

Circular and crossing movements have the same effect.

Do not hold the CryoProbe applicator still above the lesion.

It may flood the area and damage healthy tissue surrounding the lesion.

Position  
the CryoProbe correctly

This can cause blisters and pain to patients.

Condensation of humidity in the air on the frozen field of the lesion during operation may in some cases create an ice field on top of the lesion. This 'igloo of ice' will hinder the effective working of the liquid phase of the cryogen (at  $-89^{\circ}\text{C}/-128^{\circ}\text{F}$ ) and should be removed with the applicator tip.

After thawing, the tissue turns reddish and has the aspect of a zone of inflammation due to interstitial bleeding.

The learning curve for utilizing the CryoProbe is short. CryoProbe provides the possibility to have full control in the application of freezing power. Start practicing on the regular lesions and treat the more delicate ones after having acquired some experience.

### **Duration/Depth of the treatment**

Hold the tip of the applicator as close as possible to the lesion. In all literature, the rapid drop of temperature (thermal shock) is described as an essential criterion. However, always relate to the part of the body where the operation is applied. In some parts, skin is thinner which requires more caution. A typical freeze on viral infections may last from 5 to 10 seconds for a small flat wart, and up to 45 seconds for a full thickness plantar wart. For general purposes most lesions take about 2 to 30 seconds. The duration of the treatment time will depend on the surface area and the thickness of the tissue being treated. Thickness must be estimated based on clinical experience. After the first freezing cycle the tissue should be allowed to thaw for about 30 seconds followed by a second freeze. Note that the tissue will freeze faster than during the first freezing cycle.

This 'freeze-thaw-freeze' technique offers the best opportunity for success. You up visit is recommended after 2 to 4 weeks to confirm that all pathological tissue has disappeared or that a follow-up treatment is necessary.

## **Interact with your patients**

Inform patients that there will be a pain sensation similar to a ballpoint being pushed into the skin. Generally, little or no discomfort is experienced during the first few seconds when you manage to avoid treatment of healthy tissue. The ice will reach the caudal extent of the lesion after a given amount of time related to the depth of the lesion. From that moment the patient will experience a pain sensation.

This may be the moment to stop treatment.

There might be a little residual stinging for a few minutes after treatment.


Occasionally, a blister might form and persist for a few hours.

## Liability waiver

Improper use, including excess freezing levels beyond those, which are recommended or for an excessive duration, may result in bodily injury to clients/patients or to operator. H&O Equipments nv/sa and their affiliates, respective directors, officers, shareholders, employees, agents and contractors are not liable or responsible, regardless of whether such liability or responsibility is based on breach of contract, tort, strict liability, breach of warranties, failure of essential purpose, fundamental breach or otherwise, for any death or injury, whether physical or mental, or for any incidental, consequential, indirect, special or punitive damages, arising out of the CryoProbe, its design, specifications, possession and use, and treatment procedures, and whether or not any such death, injury, loss, damage(s) result from the negligence, default or error in judgment by H&O Equipments nv/sa, their affiliates, respective directors, officers, shareholders, employees, agents and contractors, and even if advised of the possibility of such damages. You agree to indemnify H&O Equipments nv/sa, their affiliates, respective directors, officers, shareholders, employees, agents and contractors, from and against any and all liability, damages, losses, costs, judgments, fines, penalties and expenses (including legal expenses) of any kind or nature, including, without limitation, incidental, consequential, indirect, special or punitive damages, arising out of claims, demands, actions, causes of action, proceeding or suits, whether in law or in equity, due to any death, injury, loss, damage or damages as hereinbefore referred to.

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